

Oral Irrigator

Installation of Water Diverter

HS0120A



1

Use the wrench to unscrew the faucet aerator body.



2

Remove the faucet aerator.



3

Choose a connector which is suitable for your faucet, and install the connector to the water diverter.



4

Be sure the washer (gasket) is on the connector, then screw it onto the end of the faucet.



5

Connect the silicone tube to the side to water diverter.

Oral Irrigator

Installation of Mouthwash Tank & Handle



6 Place a washer (gasket) between the mouthwash tank and water diverter, then connect the mouthwash bottle to the side of the water diverter.



7 Add some mouthwash, then you can rinse your teeth with mouthwash.



8 Use the Double Sided Foam Tape to stick the Combination Rack onto the wall and place the handle & the spray nozzles. Then the installation is complete.

Oral Irrigator with Mouthwash Tank - User Manual



The Switch of Water Diverter

Downward:

Tap water runs out from faucet.



The Switch of Water Diverter

Leftward:

Water sprays out from Oral Irrigator.



Bottle Up:

Rinse with mouthwash.

Bottle Down:
Rinse with tap water.

 **Reminders:**
Please put the spray nozzle into your mouth and face down before opening the faucet. Then it can prevent the water from spraying everywhere!